

Overview

	Nano & Primo	Pro	
Measurement Manoeuvres	TV	✓	
	VC	✓	
	FVC	✓	
	FIVC	✓	
	FVC + FIVC	✓	
	MVV	✓	
General Features	Post Phase (Bronchodilatation test)	✓	
	Child Incentive	✓	
	EHR Interface	✓	
	Bronchial Provocation		Optional
	Shared Database	✓	✓
	Measurement Variables	92	102
	Reference Sets	32	32
View and Report	Volume-Time Curves		✓
	Volume-Time 6s Curves		✓
	Flow-Volume Curves	✓	✓
	Best Curves	✓	✓
	Single Result		✓
	Best Result	✓	✓
	Trending	✓	✓
	Final Report	✓	✓
	Print Screen	✓	✓
	Pre/Post Report Page 1	✓	✓
	Pre/Post Report Page 2		✓
	Pre/Post Report Page 3		✓
	Pre/Post Report Page 4		✓
	PDF Export & Image Export	✓	✓



Kellolahdentie 27, FI-70460 Kuopio, Finland
 P.O. Box 54, FI-70101, Kuopio, Finland
 +358 44 777 7002
 +358 17 283 3300
 sales@medikro.com
 www.medikro.com

Your local Medikro representative: