

# MEDIKRO nano

PC-based Mobile Spirometer

- ▶ Ultimate Mobility
- ▶ Incomparable Durability
- ▶ Patient Safety
- ▶ Environmentally Friendly



Product Code: M9487

Medikro® Nano is the world's smallest Windows-compatible spirometer, weighing only 11 g. With Nano, the advantages of PC-based spirometry are now available in a mobile setting.

Nano utilizes the easy-to-use Medikro® Spirometry Software (MSS). The system features fulfil the needs of doctors who want to perform on-site testing with a mobile system, while maintaining precision and accuracy of results. With disposable Medikro® SpiroSafe flow transducers there is no need for bacterial filters or sterilization. This increases patient throughput and saves operational costs.



## Nano Highlights



### Mobility

Medikro Nano is the world's smallest Windows-compatible spirometer, weighing only 11 g. The advantages of PC-based spirometry are now available in a mobile setting.

### Small and durable

Minimized electronics reduce hardware failures, and the need for excessive maintenance procedures.



### User-friendly

Medikro Spirometry Software is user-friendly and provides intuitive navigation and use.



### Patient safety

Disposable SpiroSafe flow transducers provide security in infection control. No exhaled air gets into the pressure tube or the device. This prevents patient-to-patient cross-contamination and protects the spirometer.



### Multiple calculation criteria and interpretation options


MSS provides multiple options for calculation of Best PEF, FEV% and FEF. In addition, interpretation options are available for ventilation function, spirometry findings and the bronchodilatation test.

### PDF export

Print final reports in PDF format or as images, or let the software do it for you automatically.



## Consumables & Accessories



**MEDIKRO® SpiroSafe**

Product code: M9256



**MEDIKRO® Calibration Syringe**

Product code: M9474 (3L) / M9477 (1L)



**MEDIKRO® Ambi**

Product code: M9489



Your local Medikro representative:

🏠 Kellolahdentie 27, FI-70460 Kuopio, Finland  
✉ P.O. Box 54, FI-70101, Kuopio, Finland  
☎ +358 44 777 7002  
📞 +358 17 283 3300  
@ sales@medikro.com  
🌐 www.medikro.com