

MEDIKRO SpiroSafe

Disposable Flow Transducer

- ▶ No Filters, No Sterilization
- ▶ Patient Safety
- ▶ Cost-effective & Ecological
- ▶ Easy Twist-to-Connect System



Product Code: M9256

The Medikro® SpiroSafe flow transducer is a disposable pneumotach. With SpiroSafe there is no need for bacteria filters, sterilization or additional cardboard mouthpieces. This increases patient throughput and saves operational costs. Every single flow transducer is individually verified for accuracy.

SpiroSafe is designed to increase safety by preventing patient-to-patient cross-contamination. SpiroSafe helps to protect the environment by producing up to 78% less waste than standard bacteria filters. Medikro disposable flow transducers are made of polyethylene and can be disposed of by incineration without producing toxic emissions.



SpiroSafe Highlights



Easy Twist-to-Connect system

Medikro Twist-to-Connect makes attachment and removal of SpiroSafe flow transducers easy and fast.

No filters, no sterilization

With disposable Medikro SpiroSafe flow transducers, there is no need to use bacteria filters or to perform sterilization.



Patient safety

Disposable SpiroSafe flow transducers provide security in infection control. No exhaled air gets into the pressure tube or the device. This prevents patient-to-patient cross-contamination and protects the spirometer.



VCT technology

Vibration Control Tubing (VCT) cancels the signal noise caused by patient and tube movement. This stabilizes the signal and improves overall spirometry quality.



Cost-efficient

Medikro SpiroSafe can be changed in a matter of seconds. The whole system is then immediately clean and ready for the next patient.

Environmentally friendly

SpiroSafe helps to protect the environment by producing up to 78% less waste than standard bacteria filters. It is made of polyethylene and can be disposed of by incineration without producing toxic emissions.



Compatible With



MEDIKRO® Pro
Product code: M9488



MEDIKRO® Nano
Product code: M9487



MEDIKRO® Primo
Product code: M9492



🏠 Kellolahdentie 27, FI-70460 Kuopio, Finland
✉ P.O. Box 54, FI-70101, Kuopio, Finland
☎ +358 44 777 7002
📞 +358 17 283 3300
@ sales@medikro.com
🌐 www.medikro.com

Your local Medikro representative: