

MEDIKRO

High-Precision Volume Calibration Syringes

3000

&

1000


- ▶ High Accuracy and Precision
- ▶ Robust Aluminium Design
- ▶ Excellent Stability
- ▶ Available with Universal Adapter Kit



Product Code: M9474 (3000ml) • M9477 (1000ml)

Medikro calibration syringes are important tools for laboratory quality management. They are easy to use and meet or exceed the latest ATS/ERS 2005 recommendations for accuracy and precision.

Medikro 3000 and 1000 are made of 6061-T6 aluminium alloy to guarantee robustness and a long product lifetime. Their unique hexagonal shape, with full-length padding on each of the six sides, provides excellent stability. The syringes are designed to maintain integrity in the challenging pulmonary function testing environment.



Calibration Syringe Highlights



Accuracy & precision

Medikro calibration syringes meet or exceed all the latest ATS/ERS recommendations for accuracy and precision.

Excellent stability

The unique hexagonal design has integrated paddings on all sides. This ensures a stable calibration process and prevents the syringe from rolling.



6061-T6 aluminium alloy

Medikro calibration syringes are manufactured from the versatile 6061-T6 aluminium alloy. This guarantees robustness and ensures a long product lifetime.



Universal adapter kit

An optional universal adapter kit is available with outlet port adapters in four different sizes. This enables you to use Medikro calibration syringes with the spirometer of your choice.



Easy to use

Medikro calibration syringes are easy to use and maintain their integrity in the challenging pulmonary function testing environment.

Spirometry according to ATS/ERS

The use of calibration syringe enables calibration and verification of the spirometry system. This ensures that your system is accurate and functioning properly.



Related Products



MEDIKRO® Calibration Syringe Outlet Port Adapters

Product code: M91175-25.27.29.31

Your local Medikro representative:



🏠 Kellolahdentie 27, FI-70460 Kuopio, Finland

✉ P.O. Box 54, FI-70101, Kuopio, Finland

☎ +358 44 777 7002

📞 +358 17 283 3300

@ sales@medikro.com

🌐 www.medikro.com